


Time to pump 1+1 Liter of Glycerin + Water at constant pressure head



	30 deg	60 deg	90 deg	30 deg	60 deg	90 deg	30 deg	60 deg	90 deg	30 deg	60 deg	90 deg
	0	0	0	40 bpm	40 bpm	40 bpm	70 bpm	70 bpm	70 bpm	100 bpm	100 bpm	100 bpm
	Without	Without	Without	With Ball								
	Ball	Ball	Ball									
■ Gly+H2O at 245mm	30	10	9	32	11	9	53	14	12	53	14	12